

Handling of the Yeast Cultures

Slant Agar Cultures

Storage

- maximum 4 days at a surrounding temperature of (20 – 25°C)
- maximum 4 weeks at 4°C
- Culture medium: Wort Agar

Application

- Inoculation of the slanted agar culture in 50 ml of sterile cast wort at room temperature
- Inoculation into resp. leave onto the plant wort after 3 days in 1 l of sterile wort
- Incubation at room temperature – at the stage of rocky heads another inoculation in sterile wort in the ratio 1:10 – 1:20 at approx. 17°C
- Another inoculation in resp. leave onto the plant wort during rocky heads in the ratio 1:10 – 1:20 until an adequate pitching quantity is achieved. With every step of the propagation a gradual cooling to pitching temperature has to be done!

Liquid Culture

Storage

- maximum 1 – 2 days at a surrounding temperature of (20 – 25°C)

Application

- The bottle content (500 ml) may be inoculated directly into 10 l of sterile wort
- The increase of quantity at further steps of propagation is between 1:10 to 1:20 (for aerobe propagation)!
- With every step of the propagation a gradual cooling to pitching temperature has to be done!

For questions as to the handling of the cultures or orders please contact:

Thomas Huber, 089/85805-23; huber@doemens.org

Agnes Vels-Schmitzer, 089/85805-52; velsz-schmitzer@doemens.org